



Barkey TCP Logging Tool

Electronic Documentation Software

Documentation of date, time, temperature and the chosen program.


Records user ID and/or barcodes of samples.

Connects the Barkey plasmatherm, plasmatherm V, plasmatherm C&G and varitherm to your computer/network.

Export your data to MS Excel* and call up already saved data.

Multiple language menu.

Automatic recording of all changes done at the Barkey device (e.g. setting a new program duration).



DATA AREA

- user ID (if scanned)
- serial number of the connected Barkey plasmatherm TCP
- sample ID 1–8 (if scanned)
- IP address
- directory

PROGRAM AREA

- program view (e.g. plasma or blood)
- start, end, maximum of the recorded temperatures with date and time

EVENT AND DIAGRAM VIEW

- Live view with documentation of temperature every second (only in combination with Barkey varitherm)

► Compatible to Win 10.

► Option for IQ/OQ/PQ of TCP Logging Tool.

► Might be additional to your Quality Management

► For varitherm: temperature log filing every second.