



 **Barkey** *vialguard*

Barkey **vialguard**

Accessory for plasmatherm C&G and varitherm



▶ Silicone grid for automated dry vial thawing with Barkey plasmatherm C&G and Barkey varitherm.

▶ Conforms safely to the top and bottom warming cushions to ensure uniform heat transfer.

▶ Sleek design for tight conformity of vials helps secure positioning and prevent colliding during thawing process.

▶ Easy to write reproducible SOPs.

▶ Enables simultaneous multiple vial thawing. Supports various volumes of vials and tubes during a single thawing run.

▶ Withstands most cleanroom cleaning procedures. Compatible with common sanitizing.

Product	Barkey vialguard
Material	Silicone
Color	transparent
Dimensions	400 x 275 x 2 mm (L x W x D)
Dimensions cut out	125 x 25 mm
Weight	100 g
Certification	CE